

Respiratory Therapy

Disposable Resuscitator

Description



Adult

SW71100B Disposable Resuscitator Adult

Capacity of resuscitator	1600ml ± 10%
Material of resuscitator	PVC
Pop-off valve	60cm H ₂ O
Disposable reservoir bag	2500ml± 15%
Disposable mask	#5
Oxygen tube	2m
Air Way	80mm-100mm
Packing	Colored Box



Child

SW71200A Disposable Resuscitator Child

Capacity of resuscitator	500ml ± 10%
Material of resuscitator	PVC
Pop-off valve	40cm H ₂ O
Disposable reservoir bag	2500ml± 15%
Disposable mask	#3
Oxygen tube	2m
Air Way	60mm-80mm
Packing	Colored Box



Infant

SW71300A Disposable Resuscitator Infant

Capacity of resuscitator	280ml ± 10%
Material of resuscitator	PVC
Pop-off valve	40cm H ₂ O
Disposable reservoir bag	900ml± 15%
Disposable mask	#1
Oxygen tube	2m
Air Way	40mm-60mm
Packing	Colored Box

Respiratory Therapy

Reusable Resuscitator



Adult

SW72100A

Reusable Resuscitator Adult

Capacity of resuscitator	1600ml ± 10%
Material of resuscitator	PC, Silicon
Pop-off valve	60cm H ₂ O
Disposable reservoir bag	2500ml± 15%
Reusable mask	#3, #5
Oxygen tube	2m
Air Way	80mm-100mm
Packing	Colored Box



Child

SW72200A

Reusable Resuscitator Child

Capacity of resuscitator	500ml ± 10%
Material of resuscitator	PC, Silicon
Pop-off valve	40cm H ₂ O
Disposable reservoir bag	2500ml± 15%
Reusable mask	#2, #3
Oxygen tube	2m
Air Way	60mm-80mm
Packing	Colored Box



Infant

SW72300A

Reusable Resuscitator Infant

Capacity of resuscitator	280ml ± 10%
Material of resuscitator	PC, Silicon
Pop-off valve	40cm H ₂ O
Disposable reservoir bag	900ml± 15%
Reusable mask	#0, #1
Oxygen tube	2m
Air Way	40mm-60mm
Packing	Colored Box