

Electrosurgical Pencil Accessories

Electrosurgical Pencil Accessories

Description	Package	Ref. No.
Tip Cleaner, sterile	 =1000	SW1000003A
Tip Cleaner, non-sterile, bulk package	 =5000	SW1000003B
 *Abrasive surface effectively removes tissue *Cleans uncoated cautery blades *Adhesive backed for universal placement *X-ray detectable *Conformed RoHS requirement *single package		
Electrosurgical Holster/Quivers (HDPE), steril, 45*30*135mm	 =50	SW1000102A
		
		